



UTILITY PIPING VERIFICATION KIT

APPLICATION

To provide necessary piping components allowing utility line(s) monitoring and verification of feed water quality to STERIS products. Ensures utilities are within the recommended manufacturer specifications.

DESCRIPTION

The utility verification piping kit provides steam and water gauges, a seven function multi-meter, Technical Service and Sales Request form (water quality), and required piping components to allow monitoring on existing equipment or pre installation utility line verification.

BENEFITS

- Helps determine if performance problem is an equipment or facility issue
- Cost effective kit
- No interference to equipment operation
- Reduce unnecessary part replacement and repair cost to equipment
- Increase operational time of equipment
- Reduces install time by ensuring utilities meet operating requirements

REFERENCE LIST

Utility Piping Verification Kit	P-764332-510
<ul style="list-style-type: none"> • steam gauge 0-160 PSI • water gauge 0-100 PSI • TSSR form • miscellaneous piping parts • 10' S/S flexible hose assembly 	
Installation Guideline/Checklist	W-764332-513
Protective Case (optional)	P-764332-514

Visit **STERIS**: eservice.steris.com


NOTE: Because of STERIS's continuing program of research and development, all specifications and descriptions are subject to change without notice. Obtain certified drawings for design and installation.

For further information and order placement, please contact:

STERIS[®] STERIS Corporation
 2424 West 23rd Street
 Erie, PA 16506
 Phone: 800-333-8828
 -OR-
 Visit us online at: www.STERIS.com



Utility Piping Verification Kit
 (Typical only - some details may vary.)



Protective Case (optional)
 (Typical only - some details may vary.)