


Reduce Plastic Waste in Your Supply Chain

By offering our products in bulk delivery by specially equipped tanker trucks, STERIS helps Customers reduce plastic waste in their supply chain.

One Bulk Delivery Can Save:




12 Totes of product
or
1,752 lbs | 792 kg
of plastic



65 Drums of product
or
1,407 lbs | 641 kg
of plastic



720 Pails of product
or
2,016 lbs | 936 kg
of plastic




AT STERIS, OUR MISSION IS TO HELP OUR CUSTOMERS CREATE A HEALTHIER AND SAFER WORLD which means not only creating an environment focused on the health and safety of our people, but on the communities in which we operate. STERIS has made significant strides in reducing waste. Learn more about our sustainability initiatives. If you are interested in receiving products in bulk delivery, contact us.

Critical Process Parameters

OPTIMIZE YOUR CLEANING PROCESS



CPPs are factors to consider during each phase when developing a cleaning cycle. The dynamic between time, action, chemistry, and concentration are often referred to as “Sinner's Circle”. When one aspect of the Sinner's Circle is altered, it impacts the other parameters of the cleaning process. STERIS detergents are engineered for effectiveness on biotechnology and pharmaceutical process cleaners to improve efficiency by reducing the temperature, time and mechanical action required to clean. A cleaning process that uses a pharmaceutical detergent results in an effective, efficient and sustainable cleaning cycle.



WHY STERIS?

It is important to consider your unique system when developing a cleaning process. STERIS helps evaluate your process with unmatched technical support. STERIS offers:

An Extensive Documentation Package

To help you develop an efficient and effective cleaning process, STERIS offers an extensive documentation package for its pharmaceutical detergents such as analytical methods, toxicity studies, stability, substrate compatibility, and conductivity studies. An extensive technical library is available to support your validation objectives.

PACE™ Program

The STERIS Process and Cleaner Evaluation™ (PACE) program helps you target critical performance parameters suited for your manufacturing process, objectives, and constraints. Your process cleaner is coated on a coupon of representative material of construction. Different cleaning agents and cleaning methods are used to determine a fast and efficient cleaning recommendation.

Dedicated Scientific Support

The STERIS Technical Services team has the expertise to assist you in developing and validating a cleaning process. The STERIS Technical Services team participates in global industry events, presenting on innovative topics and trends.

Technical Sales Team

STERIS has a global direct sales team with the technical experience to assist you on-site. The STERIS sales representatives are available for cleaning trials and product support.

All company and product names are trademarks of STERIS, its affiliates or related companies, unless otherwise noted. ©2025 STERIS

