



PQ Validation Service

Performance Qualification and Requalification

Ensure your equipment meets regulatory requirements with our Performance Qualification (PQ) Validation Service. Authored and administered by our global team of subject matter experts, you gain both validation of compliance directly from the OEM and assurance that your equipment meets all of your standards, even in the most challenging applications.

Whether you're qualifying your equipment for the first time or performing periodic requalification, STERIS provides a complete solution for your applications.

PQ Validation Service includes:

- Empty Chamber Heat Distribution Study
- Heat Penetration Study
- Biological Challenge Study

Visit www.sterislifesciences.com to learn more about our complete portfolio of Validation Services.



The STERIS PQ offering is designed with you in mind. Our fundamental approach ensures that your specific application and equipment requirements are vital in developing the protocol. Together, we will scope a validation plan to meet your varying applications.

SERVICES AVAILABLE

<h3>STANDARD PQ</h3>	<ul style="list-style-type: none"> • Standard number of runs and set of tests for sterilizers and washers • Includes PQ protocol with one round of edits • On-site PQ protocol execution 	<p>Following the completion of PQ service, you receive:</p> <ul style="list-style-type: none"> • Reports and data sheets from all tests performed • Signed PQ certificate of completion • Hard copy binder containing the executed PQ protocol • Final report for your approval
<h3>CUSTOM PQ</h3>	<ul style="list-style-type: none"> • Available for sterilizers, washers, water stills, and pure steam generators • Customized to your unique applications and equipment requirements 	
<h3>REQUALIFICATION</h3>	<ul style="list-style-type: none"> • On-site testing performed according to either Customer-provided or STERIS-written protocol 	

Prior to performing PQ service, Cycle Development should be performed to determine acceptance criteria for the PQ. During Cycle Development, we help you determine load items, patterns, and cycle parameters including temperature, time, and pressure in order to meet User Specified Requirements.



STERIS CORPORATION
5960 Heisley Road
Mentor, OH 44060, USA
T +1 800 444 9009

www.STERISLifeSciences.com

STERIS Offices Worldwide

Asia Pacific	+65 (0) 68 41 7677
Canada	+1 800 444 9009
Finland	+358 9 25851
France	+33 (0) 2 38 84 85 40
Germany	+49 (0) 221 466120
Italy	+39 (0) 2 21303 424
Latin America	+1 800 444 9009
Spain	+34 (0) 916 585 920
United Kingdom	+44 (0) 1256 840400

STERIS Life Sciences

STERIS Life Sciences is your trusted partner in contamination control.

For over 100 years, STERIS has been a global leader and expert in the industries of sterilization, cleaning and infection control. Today, we continue to build on this heritage by providing the products and resources you need throughout every step of your process. From innovation to results, our comprehensive portfolio of products and services is designed to meet your needs.

With offices located worldwide and technicians in over 100 countries, you can count on us to be where you are when you need us, every step of the way. Trust STERIS Life Sciences to help you create a healthier and safer world.

© 2018 STERIS Corporation
All rights reserved

No part of this publication may be reproduced in any material form without the permission of STERIS Corp.