

# Services to Consider When Validating Your Equipment

Complying with regulatory requirements is critical to your success. Whether you're qualifying new equipment or performing regular requalification, having a robust validation plan provides you with documentation and assurance that your equipment operates up to your standards.

Validation services can help prevent costly downtime and keep you in compliance. With procedures and protocols created based on current Good Manufacturing Practices (cGMPs) and backed by over a century of technical expertise, STERIS offers comprehensive validation packages throughout your equipment's lifecycle.

1

### Factory Acceptance Testing (FAT)

Verify your equipment has been built to your specific design qualifications with comprehensive testing in our manufacturing facilities.

2

#### Site Acceptance Testing (SAT)

Ensure your equipment works to your specifications in your facility. Support your qualification process with comprehensive testing and documentation.

3

# Installation Qualification (IQ) & Operational Qualification (OQ)

Confirm your equipment was installed correctly and operates according to your specifications with equipment-specific documents and on-site execution.

4

# **Cycle Development**

Enhance your
equipment's performance
by defining the optimal
load configuration and
customizing cycle
parameters for your
equipment.

5

## Performance Qualification (PQ) & Requalification

Validate that your equipment consistently meets your performance specifications with documented protocols and testing.

